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GB Instructions for Use

Always carefully read the instructions for use before using NOBASPONGE®-border, as the information is updated on a regular basis. The product may be used by laypersons under supervision of healthcare professionals.

Product description:

NOBASPONGE®-border is a sterile, absorbent foam dressing featuring a hydrocellular pad with a low-friction, breathable polyurethane film backing and an atraumatic silicone wound contact layer. The product can be used as a low-trauma primary dressing on fragile skin to absorb wound exudate and provide an optimal wound healing environment, even under compression.

Components:

Polyurethane film, silicone gel, polyurethane foam, release film

Indication:

- Chronic and acute wounds with light to moderate levels of exudate
 - Skin abrasions
 - Pre-tibial lacerations
 - Pressure ulcers
 - Diabetic foot ulcers
 - Leg ulcers
 - Post-operative wounds
 - Superficial burns
 - Donor sites
 - Traumatic wounds
- Shallow, granulating wounds

Contra-indications:

Arterial bleeds, allergies against the material

Warnings and precautions:

- NOBASPONGE®-border is not suitable for arterial bleeds.
- Do not overlay multiple foam dressings or overlap them
- Do not stretch the skin or dressings during application. Silicone is hydrophobic and will not adhere to wet skin.
- Do not use with oxidizing solutions such as hypochlorite or hydrogen peroxide.
- Before using a sterile product, visually inspect the packaging to ensure that it is intact. Do not use if individual package is damaged or open.
- Single use device: The reuse of single use devices may lead to a microbiological hazard. Treatment for reuse may adversely affect the performance of the product.

Should serious incidents occur in relation to the device, they are, according to MDR (EU) 2017/745, to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Application instructions:

- Cleanse the wound in accordance with local clinical protocols.
- Select an appropriate-sized dressing. Ensure the foam pad covers the entire wound area and covers a minimum of 2 cm undamaged skin beyond the wound edges.
- Prepare and cleanse the skin surrounding the wound area and remove excess moisture.
- Remove NOBASPONGE®-border from the packaging using a clean technique.
- Remove the release liners and apply pad-side down to the wound.
- Smooth the dressing border down to ensure good adherence with the peri-wound skin.
- Always ensure the dressing is not applied too tightly.

Frequency of Dressing Change:

- The use of the dressings should always be re-evaluated if wound exudate levels change.
- During the early stages of wound management, observe the wound daily for signs of infection (it can be signalled by the occurrence of higher body temperature, pain, redness, swelling around the area) or other complications (or as recommended by local clinical protocols).
- Infected wounds should be inspected and treated as per local clinical protocol.
- The frequency of dressing change is determined by the condition of the wound and local clinical protocols.
- The dressing can be left in place, undisturbed, for up to 7 days, or until exudate is visible and approaches 2 cm from the edge of the dressing, whichever is sooner. When the amounts of exudates become smaller or when there are no exudates, the dressing should be changed as often as needed so that it does not adhere to the wound.

Dressing Removal:

Careful and correct dressing removal is vital to minimize the risk of damaging the wound bed and new granulating tissue; please wear gloves.

- Remove the retention dressing or bandage and carefully lift the dressing away from the wound.
- The dressing should be removed slowly and carefully following the direction of hair growth.
- As the dressing is removed, the newly exposed skin should be supported.

Side effects:

Unknown

Storage and transport:

The product has a shelf-life of 3 years. The product must no longer be used after the expiry date has passed. The product should be stored horizontally.

Waste disposal:

In accordance with the requirements of the applicable local legislation and with the standards for the prevention of infections.

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GB NOBAMED Paul Danz AG