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<span><span><span></span></span></span> EN <span> </span> INSTRUCTION FOR USE - PEG Dressing Kit
<b>Product name:</b> Dressing Kit
<b>Material:</b> <ul style="list-style-type: none"><li>Component: 2x Non-woven swabs 7.5x7.5cm 4-ply, 2 x Non-woven split swabs 7.5x7.5cm 4-ply, 1 x Elastic plaster (round corners) 12x15cm.</li></ul>
- Elastic plaster is made from PE, non-woven is made from viscose and polyester.
- Sterile packed, quick and easy to use.

**Intended purpose:**
Dressing kit is used for percutaneous endoscopic gastrostomy (PEG)/percutaneous endoscopic jejunostomy (PEJ).

**Intended users:**
The product shall be used by or under the supervision of a healthcare professional.

**The combination of the device with product:**
Not applicable

**Contraindications:**
People who are allergic to product materials

**Complications or side-effects:**
Unknown

**Instruction of use:**

- Create a clean working surface
- Hygienic hand disinfection, put on hygienic gloves
- Remove old dressing ( in case of probe care open fixation of the holding plate and pull backwards)
- Change gloves
- Cleans skin with sterile non-woven compress disinfecting
- Clean holding plate, probe with the second compress disinfecting
- Allow to dry disinfecting area
- Put sterile non woven slit compress around the probe. If available position between skin and holding plate
- Cover holding plate with third non-woven
- Fix the dressing with the fixation plaster without tension on the skin

- Attention:**
- For single-use only.
  - Do not use if the packaging is damaged.
  - Mobilisation can of the tube(PEG/PEJ) should be carried out according to the tube manufacturer’s recommendations.
  - Mobilisation can prevent the inner retaining plate from growing into the inner wall of the stomach (buried bumper syndrome).
  - Do not use on the patients with a known hypersensitivity to the product itself or to its components.
    - Stop using immediately in case of allergy.
  - It will be treated as medical waste after used, Please refer to the local/national regulations regarding waste disposal.

**Storage:**
Store in rooms with a stable temperature from 0 °C - 45 °C and humidity from 10% - 70%.

<span><span></span></span>	Medizinprodukt / Medical device / Dispositif médical / Dispositivo médico / Dispositivo medico / Medisch hulpmiddel / Medicinsk udstyr
<span><span></span></span>	Eindeutige Gerätekennzeichnung / Unique device identifier / Identifiant unique du dispositif / Identificador único del dispositivo / Identificatore unico del dispositivo / Uniek apparaatnummer / Unik enhedsidentifikator
<span><span></span></span>	Sterilisiert mit Ethylenoxid / Sterilized using ethylene oxide / Stérilisé avec de l’oxyde d’éthylène / Esterilizado utilizando óxido de etileno / Sterilizzato con ossido di etilene / Gesteriliseerd met behulp van ethyleenoxide / Steriliseret med ethylenoxid
<span><span></span></span>	Nicht erneut sterilisieren / Do not resterilize / Ne pas restériliser / No esterilizar / Non risterilizzare / Niet opnieuw steriliseren / Skal ikke steriliseres igen
<span><span></span></span>	Einfach-Steril-Barrieresystem / Single sterile barrier system / Système de barrière stérile / Sistema de barrera estéril único / Sistema di barriera sterile semplice / Enkelvoudig steriel barriersysteem / Enkelt sterilt barriersystem
<span><span></span></span>	Einfach-Steril-Barrieresystem mit äußerer Schutzverpackung / Single sterile barrier system with protective packaging outside / Emballage de protection non stérile avec, à l’intérieur, un système de barrière stérile / Sistema de barrera estéril único con un embalaje de protección exterior / Sistema singolo di barriera sterile con imballaggio protettivo esterno / Enkelvoudig steriel barriersysteem met beschermende buitenverpakking / Enkelt sterilt barriersystem med indre beskyttelses forpakning

Keep away from direct sunlight and keep dry.
Store in a non-corrosive gas environment, away from light, moisture, fire and flammable materials.

**Lot no.:** see labelling.

**Date of manufacture:** see labelling.

**End of manufacture:** see labelling.

**Validity:** 3 years.

**Device related adverse event reporting:**
Any adverse event involving the Dressing Kit manufactured by Kingphar should be reported to Kingphar medical or local European representative immediately.

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Hubei Qianjiang Kingphar Medical Material Co., Ltd. Yuanguang Road   433100 Qianjiang   PEOPLE'S REPUBLIC OF CHINA
Shanghai International Holding Corp. GmbH (Europe)   Eiffestrasse 80   20537 Hamburg   Germany
CMC Medical Devices GmbH   Bahnhofstrasse 32   CH-6300 Zug   Switzerland