

# SPECIFICATIONS

## ROMED ANEROID SPHYGMOMANOMETERS

### REF: BPM-AN-50



*Applicable harmonized standard: EN ISO 81060-1:2012*  
*Medical device conform MDD 93/42/EEC*

Aneroid sphygmomanometer with adult sized cuff, with separate meter and bulb. A diagnostic apparatus for taking human blood pressure. Small in size and easy to carry. The diagnostic efficiency of the apparatus is similar to that of the mercurial type.

#### **Specifications**

Measuring range : 0-300 mmHg  
Accuracy :  $\pm 2.0$  mmHg

Cuff size : 50 x 14 cm, one-size-fits-all for adults  
Cuff material : Nylon with internal latex bladder and Velcro closure system, washable  
Cuff colour : Blue, with artery indicator and gauge holder

Inflation bulb : Manual, with manual air release valve  
Tubes length : 50 cm  
Bulb, tube material : Rubber, black

Meter : Gauge to allow recalibration

#### **Precautions**

- When using this apparatus, do not let the pointer exceed 300 mmHg by inflation, because if the internal spring exceeds the elastic limit, accuracy can be affected.
- Avoid violent shocks in order to keep the reading accurate.
- Put the meter in the black vinyl carrying case after use and always keep it clean.

#### **Packing**

Per piece in an inner box, 50 pieces in an export carton.